
PakGent[®]

PRODUCT INTRODUCTION

Centrifuge Tubes (15ml & 50ml)

PakGent Bioscience (Suzhou) Co., Ltd.

Tel: +86 (0) 512 6670 0017


Email: info@pakgentbio.com

Centrifuge tubes are widely used in laboratories for a variety of purposes, including but not limited to storing samples and performing high-speed centrifugation.

- USP Class VI raw material (PP)
- SAL=10⁻⁶
- Centrifugal force: 15000 xg (conical bottom)
3500 xg (self-standing)
- rack pack & bulk pack are available.



01 Threaded cap for one hand operation.

02 User-friendly with a large writing area. Smooth writing.

03 Uniform wall thickness and stable centrifugal effect.

04 Ergonomic, collapsible storage racks made from recycled materials are better for the lab and the environment. Available for both 15 ml and 50 ml sizes.

DECLARATION OF COMPLIANCE

DOC No. PK20240204001 Issued Date: 4th, February, 2024

This certificate of conformity is based on an evaluation of a sample of the following mentioned product, and technical report and documentation are at the license holder's disposal. We herewith declare that the above mentioned products meet the provisions of the Council Directive 93/42/EEC concerning Medical Device, as amended by 2007/147/EC. This declaration of conformity is issued under the sole responsibility of the manufacturer.

Manufacturer: PakGent Bioscience (Suzhou) Co., Ltd.
Manufacturing site: Floor 1-3, Building 33, No. 226, Gongtang Road, Luzhi Town, Wuzhong District, Suzhou, China (215127)
Brand: PakGent
Product: Centrifuge Tubes
Par.No: CT50R, CT50B, CT50RS, CT50BS, CT15R, CT15B, MCT50, MCT20, MCT15, MCT05


Note: The Declaration of Conformity, all document and certificates and attestations shall be duly certified true copy by the Applicant.

Authorised Signatory:
 授权签字人: Alvin Qag
 Position: QA Manager
 Date: 4th, February, 2024

PakGent PakgentBio(Changshu)Co.Ltd

9. Test results
 9.1 Observing images under a microscope

9.2 Result determination

Group	Replicate1	Replicate2	Replicate3	Reactivity	Cytotoxicity
Sample	0	0	0	0	None
Blank Control	0	0	0	0	None
Negative Control	0	0	0	0	None
Positive Control	4	4	4	4	Severe

10. Conclusion
 Under the conditions of this experiment, the cytotoxicity test results of the sample are qualified.

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
Page 3 of 4
 PROPRIETARY CONFIDENTIAL - INTERNAL USE ONLY
 仅供内部使用, 请勿外传


Brand	PakGent		Global Brands						Local Brands					
			F-		N-		G-		J-		M-		C-	
Model	15ml	50ml	15ml	50ml	15ml	50ml	15ml	50ml	15ml	50ml	15ml	50ml	15ml	50ml
Centrifugal force (g)	15,000	15,000	12,000	16,000	11,000	17,000	15,000	17,000	13,000	12,000	13,000	13,000	13,000	12,000
Leak detection 70kpa for 2min	OK	OK	NG	NG	OK	OK	OK	OK	NG	OK	OK	NG	OK	NG
Touch of screw cap after high temperature and high pressure	9	8	9	5	8	9	8	8	8	7	8	6	7	6
<p>Note: Feel score:</p> <p>9-10 points can easily screw the lid, feel silky.</p> <p>7-8 points can be normal capping, slightly dry.</p> <p>6 points can screw the cap, but it is quite tight.</p> <p>5 points No capping.</p>														

Ordering Information


Part No.	Volume	Centrifuge Force	Bottom	Package Type	Pack
CT50B	50ml	15000 xg	conical bottom	bulk pack	25 pcs/bag, 20 bags/case
CT50R	50ml	15000 xg	conical bottom	rack pack	25 pcs/rack, 20 racks/case
CT15B	15ml	15000 xg	conical bottom	bulk pack	25 pcs/bag, 20 bags/case
CT15R	15ml	15000 xg	conical bottom	rack pack	25 pcs/rack, 20 racks/case
CT50BS	50ml	3500 xg	self-standing	bulk pack	25 pcs/bag, 20 bags/case
CT50RS	50ml	3500 xg	self-standing	rack pack	25 pcs/rack, 20 racks/case

15ml conical bottom


50ml conical bottom

50ml self-standing

PakGent®

THANKS!

PakGent Bioscience (Suzhou) Co., Ltd.

Tel: +86 (0) 512 6670 0017

Email: info@pakgentbio.com