

Centrifuge Tubes

PakGent®


Guidelines For Use

No.: PGCT20240401E




PakGent® Centrifuge tubes are designed for your everyday centrifugation and storage needs. These tubes include large white writing area and black easy to read graduations. Racked tubes are packaged in a reusable, autoclavable, and freezable (-80°C) polypropylene rack that inter-locks to allow you to create a workstation.

Features	Description
Material	Polypropylene (PP), premium surface
Centrifugal force	Centrifuge Tubes: 15000 xg, microcentrifuge tubes: 25000 xg
Temperature range	Autoclavable (121 °C) and freezable (-80 °C)



Uniform wall thickness and stable centrifugal effect.




Self-sealing bag design, more convenient for laboratory use.



Ergonomic, collapsible storage racks made from recycled materials are better for the lab and the environment. Available for both 15 mL and 50 mL sizes.



Part No.	Volume	Centrifuge Force	Bottom	Package Type	Pack
CT50B	50ml	15000 xg	conical bottom	bulk pack	25 pcs/bag, 20 bags/case
CT50R	50ml	15000 xg	conical bottom	rack pack	25 pcs/rack, 20 racks/case
CT15B	15ml	15000 xg	conical bottom	bulk pack	25 pcs/bag, 20 bags/case
CT15R	15ml	15000 xg	conical bottom	rack pack	25 pcs/rack, 20 racks/case
CT50BS	50ml	3500 xg	self-standing	bulk pack	25 pcs/bag, 20 bags/case
CT50RS	50ml	3500 xg	self-standing	rack pack	25 pcs/rack, 20 racks/case



50ml conical bottom

50ml self-standing


15ml conical bottom

Microcentrifuge tubes




Microcentrifuge tubes are made from an ultra-clear polypropylene resin that allows for easy viewing of contents and fit most common rotors for near-universal usage. Caps can either be screw-on or attached to the microcentrifuge tubes for easy popping on, and even pop-on caps are immune to popping off during intense centrifugation.


Part No.	Volume	Centrifuge Force	Package Type	Pack
MCT50	5.0ml	25000 xg	bulk pack	200 pcs/bag, 10 bags/case
MCT20	2.0ml	25000 xg	bulk pack	500 pcs/bag, 10 bags/case
MCT15	1.5ml	25000 xg	bulk pack	500 pcs/bag, 10 bags/case
MCT05	0.5ml	25000 xg	bulk pack	1000 pcs/bag, 10 bags/case




MCT50



MCT20



MCT15



MCT05

PakGent[®]

PakGent Bioscience (Suzhou) Co., Ltd.

Address: No. 9 Lecheng Road, Bixi Street, Changshu, Jiangsu Province, China

Tel: +86 (0) 512 6670 0017

Email: info@pakgentbio.com

Website: www.pakgentbio.com



Catalog